

2018

Guidebook Endoscopic Systems

including accessories



FUJIFILM
Value from Innovation

ELUXEO™ 700 SERIES ENDOSCOPES

	Model	Observation range	Field of view	Ø Distal end	Ø Insertion tube	Ø Working channel	Bending capability	Bending capability	Working length	Image area, forceps & water jet entry	Remarks
		(mm)		(mm)	(mm)	(mm)	(up/down)	(left/right)	(mm)		
 	ELUXEO® EG-760R	2–100	140°	9.2	9.3	2.8	210°/90°	100°/100°	1,100		LCI, BLI Water jet
	ELUXEO® EG-760Z	Normal 3–100 Close 1.5–2.5	Normal 140° Close 56°	9.9	9.8	2.8	210°/90°	100°/100°	1,100		LCI, BLI Multi Zoom*1 Water jet
	ELUXEO® EG-760CT	2–100	140°	10.5	10.8	3.8	210°/90°	100°/100°	1,100		LCI, BLI Water jet
	ELUXEO® EG-740N	3–100	140°	5.8	5.9	2.4	210°/90°	100°/100°	1,100		LCI, BLI

Colonoscopes

 	ELUXEO® EC-760R-VM/VI/VL	2–100	170°	12.0	12.0	3.8	180°/180°	160°/160°	1,330/1,520/ 1,690		LCI, BLI Water jet "ColoAssist Adjust"
	ELUXEO® EC-760ZP-VM/VL	Normal 3–100 Close 1.5–2.5	Normal 140° Close 56°	11.7	11.8	3.2	180°/180°	160°/160°	1,330/1,690		LCI, BLI Multi Zoom*1 Water jet "ColoAssist Adjust"
	ELUXEO® EC-740TM/TL	3–100	140°	9.8	10.7	3.2	210°/160°	160°/160°	1,330/1,990		LCI, BLI Water jet

 To learn more about LCI and BLI visit www.bli.eu

600 SERIES ENDOSCOPES

	Model	Observation range	Field of view	Ø Distal end	Ø Insertion tube	Ø Working channel	Bending capability	Bending capability	Working length	Image area, forceps & water jet entry	Remarks
		(mm)		(mm)	(mm)	(mm)	(up/down)	(left/right)	(mm)		
	EG-600WR	2–100	140°	9.2	9.3	2.8	210°/90°	100°/100°	1,100		Water jet

Colonoscopes

	EC-600WM/WI/WL	2–100	140°	12.0	12.0	3.8	180°/180°	160°/160°	1,330/1,520/ 1,690		Water jet
--	----------------	-------	------	------	------	-----	-----------	-----------	-----------------------	---	-----------

*1 Multi Zoom: Optical magnification x 135.

For information on the product availability in your area, please contact your local partner.

500 SERIES ENDOSCOPES

Model	Observation range	Field of view	Ø Distal end	Ø Insertion tube	Ø Working channel	Bending capability	Bending capability	Working length	Image area, forceps & water jet entry	Remarks
	(mm)		(mm)	(mm)	(mm)	(up/down)	(left/right)	(mm)		

Ultra-slim gastroscopes

	EG-580NW2	3–100	140°	5.8	5.9	2.4	210°/90°	100°/100°	1,100		
	EG-530NP	3–100	120°	4.9	5.1	2.0	210°/120°	–/–	1,100		

Gastroscopes

	EG-580RD	3–100	140°	9.8	9.8	3.2	210°/120°	100°/100°	1,100		Water jet
	EG-530WR	4–100	140°	9.4	9.3	2.8	210°/90°	100°/100°	1,100		
	EG-530CT	3–100	140°	10.8	10.8	3.8	210°/90°	100°/100°	1,100		Water jet
	EG-530D	3–100	140°	11.5	11.5	2.8/3.8	210°/90°	100°/100°	1,090		Dual channel Water jet

Duodenoscopes

	ED-530XT	4–60	100° (retro 8°)	13.1	11.5	4.2	130°/90°	90°/110°	1,250		
	ED-530XT8	4–60	100° (retro 8°)	13.1	11.5	4.2	130°/90°	90°/110°	1,250		

Colonoscopes

	EC-580RD M/L	3–100	140°	9.8	10.5	3.2	210°/160°	160°/160°	1,330/1,690		Water jet
	EC-530WM3/WI3/WL3	3–100	140°	12.8	12.8	3.8	180°/180°	160°/160°	1,330/1,520/ 1,690		Water jet
	EC-530MP	3–100	140°	11.0	11.1	3.2	180°/180°	160°/160°	1,330		
	EC-530LP	3–100	140°	11.0	11.1	3.2	180°/180°	160°/160°	1,690		
	EC-530DM/DL	3–100	140°	12.8	12.8	3.8/2.8	180°/180°	160°/160°	1,330/1,690		Dual channel Water jet
	ES-530WE	3–100	140°	12.8	12.8	3.8	180°/180°	160°/160°	790		Water jet



500 SERIES ENDOSCOPES

Model	Observation range	Field of view	Ø Distal end	Ø Insertion tube	Ø Working channel	Bending capability	Bending capability	Working length	Image area, forceps & water jet entry	Remarks
	(mm)		(mm)	(mm)	(mm)	(up/down)	(left/right)			

Double-Balloon Endoscopes

	EN-580T	2–100	140°	9.4	9.3	3.2	180°/180°	160°/160°	2,000		*2
	EN-580XP	2–100	140°	7.5	7.7	2.2	180°/180°	160°/160°	2,000		*3

Short Double-Balloon Endoscope

	EI-580BT	2–100	140°	9.4	9.3	3.2	180°/180°	160°/160°	1,550		*4
--	----------	-------	------	-----	-----	-----	-----------	-----------	-------	---	----

Please check our "Guidebook DBE Accessories" for details on accessories

Bronchoscopes

	EB-580S	2–100	120°	5.3	5.1	2.2	210°/130°	–/–	600		*5 LCI, BLI HF-compatible LASER blocking filters (*6)
	EB-580T	2–100	120°	5.8	5.9	2.8	180°/130°	–/–	600		*5 HF-compatible LASER blocking filters (*7)
	EB-530H	3–100	140°	5.4	4.9	2.0	180°/130°	–/–	600		*5 HF-compatible LASER blocking filters (*6)
	EB-530S	3–100	120°	4.9	4.9	2.0	180°/130°	–/–	600		*5 HF-compatible LASER blocking filters (*6)
	EB-530T	3–100	120°	5.8	5.9	2.8	180°/130°	–/–	600		*5 HF-compatible LASER blocking filters (*6)
	EB-530P	3–100	120°	3.8	3.8	1.2	180°/130°	–/–	600		*5 LASER blocking filters (*6)
	EB-530XT	3–100	120°	6.2	6.3	3.2	180°/130°	–/–	600		*5 HF-compatible LASER blocking filters (*6)

*2 Balloon controller PB-30, balloon BS-4 and overtube TS-1314B are required. Both products are latex-free.

*3 Balloon controller PB-30, balloon BS-4 and overtube TS-1114B are required. Both products are latex-free.

*4 Balloon controller PB-30, balloon BS-2 and overtube TS-13101 are required.

*5 Suction button SB-602 is required.

*6 Diode LASER (810 nm).

*7 Nd-Yag LASER (1064 nm) & Diode LASER (810 nm).



ENDOSCOPY SYSTEMS

Model	Product Description
-------	---------------------

Video Processors and Light Sources

ELUXEO™ VP-7000 BL-7000 HD Full HD endoscopy		▶ 4-LED Multi Light™ system ▶ BLI and LCI modes ▶ Full HD 1920 x 1080p digital output ▶ DICOM compatible ▶ Applicable to USB Flash memory ▶ DK-7000E keyboard ▶ Low energy, long-lasting LED light source ▶ Image output connectors: HD-SDI x 2, DVI-D x 2, DVI-I, RGB TV, S-Video, Video	Video processor with light source for 500, 600 and 700 series endoscopes ▶ VP-7000 processor 9.0 kg, 390 mm (W) x 110 mm (H) x 485 mm (D) ▶ BL-7000 light source 12 kg, 390 mm (W) x 155 mm (H) x 485 mm (D)
--	---	---	---

 To learn more about LCI and BLI visit www.bli.eu

EPX-3500HD HD Full HD endoscopy		▶ 300 W Xenon lamp ▶ Full HD 1920 x 1080p digital output ▶ DICOM compatible ▶ Applicable to USB Flash memory ▶ DK-4450E keyboard ▶ 3 patterns FICE installed ▶ Image output connectors: DVI x 2, RGB TV, S-Video, Video	Video processor with light source for 500 and 600 series endoscopes ▶ EPX-3500HD: - VP-3500HD processor 8.0 kg, 390 mm (W) x 105 mm (H) x 460 mm (D) - XL-4450 light source 15 kg, 390 mm (W) x 155 mm (H) x 450 mm (D)
--	---	---	---

*1 Ultrasonic videoscopes EG-530UR and UT cannot be connected



EUS SYSTEMS

Model	Observation range	Field of view	Ø Distal end	Ø Insertion tube	Ø Working channel	Bending capability	Bending capability	Working length	Frequency	Scanning method	Remarks
	(mm)		(mm)	(mm)	(mm)	(up/down)	(left/right)	(mm)	(MHz)		

Ultrasonic videoscopes

	EG-580UR ¹	3–100	140°	11.4	11.5	2,8	190°/90°	100°/100°	1,250	5, 7.5 10, 12	Electronic radial 360°	SU-1
	EG-580UT ²	3–100	140° (forward / oblique 40°)	13.9	12.4	3.8	150°/150°	120°/120°	1,250	5, 7.5 10, 12	Electronic curved linear array 150°	SU-1
	EB-530US ³	3–100	120° (forward / oblique 10°)	6.7	6.3	2.0	130°/90°	–/–	610	5, 7.5 10, 12	Electronic Curved linear array 65°	SU-1

Type	Model	Working length	Ø Outer diameter	Frequency	Suitable working channel	Scanning method	Remarks
		(mm)	(mm)	(MHz)			

Ultrasonic mini probes

P-series for Bronchoscopy	PB2020-M	2,150	1.4 (distal) – 1.9 (proximal)	20.0	Larger than 2.0 mm	Mechanical radial	SP-900
P-series for Gastroenterology	P2620-M	2,200	2.6	20.0	Larger than 2.8 mm	Mechanical radial	SP-900
	P2615-M	2,200	2.6	15.0	Larger than 2.8 mm	Mechanical radial	SP-900
	P2612-M	2,200	2.6	12.0	Larger than 2.8 mm	Mechanical radial	SP-900
	P2620-L	2,700	2.6	20.0	Larger than 2.8 mm	Mechanical radial	SP-900
	P2615-L	2,700	2.6	15.0	Larger than 2.8 mm	Mechanical radial	SP-900
	P2612-L	2,700	2.6	12.0	Larger than 2.8 mm	Mechanical radial	SP-900

Model	Product Description
-------	---------------------

Ultrasonic processor

SU-1		▶ 13.0 kg, 390 mm (W) x 135 mm (H) x 485 mm (D), AC 100 – 240 V ▶ Flat CP-1 keyboard with touch pad or trackball ▶ Probe type: Radial, Curved linear array ▶ Ultrasound scanning mode: 5–12MHz/Power Doppler / Colour Doppler / PW Doppler / M mode / B mode / TH/CH/CHI/Elastography ▶ Available as SU-1-S- and SU-1-H-. CHI and Elastography modes are available only in SU-1-H-
------	---	--

Ultrasonic mini probe processor

SP-900		▶ 8 kg, 377 mm (W) x 80 mm (H) x 480 mm (D), AC 100 – 240 V ▶ CP-900 control pad ▶ Scanning method: mechanical radial ▶ Ultrasound scanning mode: B mode
--------	---	--

¹ B20UR balloon is necessary
² B20UT balloon is necessary
³ B20BU balloon is necessary



FIBRESCOPIES

Model	Observation range	Field of view	Ø Distal end	Ø Insertion tube	Ø Working channel	Bending capability	Bending capability	Working length	Size of image	Remarks
	(mm)		(mm)	(mm)	(mm)	(up/down)	(left/right)	(mm)		
Gastrosopes										
FG-1Z	3–100	105°	9.5	9.8	2.8	210°/90°	100°/100°	1,030		
FG-1ZP	3–100	105°	7.8	7.9	2.2	210°/90°	100°/100°	1,030		
Colonoscopes										
FC-1Z	3–100	120°	13.0	13.0	3.2	180°/180°	160°/160°	1,660		
Bronchoscopes										
FB-120S	1–50	120°	4.8	4.9	2.2	180°/130°	–/–	600		HF-compatible
FB-120T	1–50	120°	5.9	5.9	2.8	180°/130°	–/–	600		HF-compatible
FB-120P	1–50	100°	2.7	2.8	1.2	180°/130°	–/–	600		
FB-120MP	1–50	120°	4.8	4.9	2.2	180°/130°	–/–	600		Mobile type + HF-compatible

Model	Product Description
-------	---------------------

Light sources

PS2-HP	 Light source for all fibrescopes ▶ Air pump ▶ 150 W halogen lamp ▶ 8.6 kg ▶ 195 mm (W) x 215 mm (H) x 350 mm (D)
--------	--

SPECIAL ENDOSCOPES

Model	Observation range	Field of view	Ø Distal end	Ø Insertion tube	Ø Working channel	Bending capability	Bending capability	Working length	Size of image	Remarks
	(mm)		(mm)	(mm)	(mm)	(up/down)	(left/right)	(mm)		
ENT/Naso-Pharyngo-Laryngoscopes										
Fibre	FR-120F	3–50	90°	3.3	3.2	–	130°/130°	–/–	300	
Video	ER-530S2	3–50	90°	2.9	3.2	–	130°/130°	–/–	300	110%
Video	ER-530T	2–50	140°	5.4	4.9	2.0	130°/130°	–/–	400	130%



ESD ACCESSORIES

Model	Identifier	Knife length (mm)	Working length (mm)	Ø Suitable working channel	Remarks
FlushKnife					
NEW FlushKnife CTS	DK2620J	- B15S -	2,000	2.8	Ball tip Slim Type
		- B20S -			
		- B25S -			
		- B30S -			
NEW FlushKnife INS	DK2620J	- N10S -	2,000	2.8	Needle tip Slim Type
		- N15S -			
		- N20S -			
		- N25S -			
		- N30S -			
FlushKnife BT	DK2618J	- B15 -	1,800	2.8	Ball tip
		- B20 -			
		- B25 -			
		- B30 -			
	DK2623J	- B15 -	2,300	2.8	Ball tip, Long type
		- B20 -			
FlushKnife	DK2618J	- N10 -	1,800	2.8	Needle tip
		- N15 -			
		- N20 -			
		- N25 -			
		- N30 -			
	DK2623J	- N15 -	2,300	2.8	Needle tip, Long type
		- N20 -			

Model	Identifier	Knife length	Working length	Ø Suitable working channel	Remarks
		(mm)	(mm)		
ClutchCutter					
	DP2618DT - 35 -	- 35 -	3.5	1,800	2.8
	DP2618DT - 50 -	- 50 -	5.0	1,800	2.8



ESD ACCESSORIES

Model	Outer diameter (mm)	Inner diameter of tip (mm)	Tip length (mm)	Accessory guide	Applicable endoscope	Remarks
ST hood for ESD and POEM						
DH-15GR	12.2	7.0	8.3	Yes	EG-760R, EG-760Z, EG-760CT, EC-740T, EG-600WR, EG-600ZW, EG-590WR, EG-530WR, EG-580RD	
DH-16CR	16.0	7.0	8.3	Yes	EC-760ZP-VM/VL, EC-760R-VM/VI/VL, EC-600WM, WI, WL, EC-600ZW M/L, EC-590WM4, WI4, WL4, EC-590ZW3 M/L, EC-530WM3, WI3, WL3	
DH-28GR	11.8	8.0	7.0	Without	EG-760Z, EG-760R, EG-600WR, EG-600ZW, EG-590WR, EG-530WR, EG-580RD	Short type
DH-29CR	13.0	8.0	7.0	Without	EG-760CT, EC-740T, EG-590ZW, EC-580RD, EC-530MP, EC-530LP	Short type
DH-30CR	14.8	8.0	7.0	Without	EC-760R-VM/VI/VL, EC-760ZP-VM/VL, EC-600WM, WI, WL, EC-600ZW M/L, EC-590WM4, WI4, WL4, EC-590ZW3 M/L, EC-530WM3, WI3, WL3	Short type

Observation hood

DH-11GZ	13.2		2.5		EG-590ZW	
DH-13CZ	15.7		3.0		EC-760ZP-VM/VL, EC-760R-VM/VI/VL, EC-600WM, WI, WL, EC-600ZW M/L, EC-590WM4, WI4, WL4, EC-590ZW3 M/L, EC-530WM3, WI3, WL3	
DH-14EN	10.8		1.5		EG-760Z, EG-600ZW	
DH-17EN	11.5		1.5		EI-580BT, EN-580T	

Model	Applicable endoscope	Remarks
Top Endoscopic Hoods*		
SHM	EG-760Z	Transparent type 5 pcs / box
SHL	EC-760ZP	Transparent type 5 pcs / box
SHM-B	EG-760Z	Black type 5 pcs / box
SHL-B	EC-760ZP	Black type 5 pcs / box



New submucosal injection agent designed for easy and safe resection procedures.



*Manufactured by Top Corporation, Japan.
Not for sale in Canada.

Manufactured for Cosmo Technologies Ltd., Ireland.
Not for sale in Canada.

PERIPHERAL DEVICES

Model		Product Description
PB-30		Balloon controller for DBE endoscopes ▶ 7 kg (main unit), 145 mm (W) x 170 mm (H) x 410 mm (D), AC 100-230 V ▶ Pressure of balloon: 5.6 kpa ± 2 kpa ▶ Max. flow rate of pump: 170 ± 50 ml/10 sec.
JW-2		Water pump for FUJIFILM endoscopes with a water jet channel ▶ Foot switch operation ▶ Flow late adjustment function ▶ 2.6 kg ▶ 145 mm (W) x 130 mm (H) x 265 mm (D): Main unit
GW-100		CO2 gas regulator for FUJIFILM endoscopes ▶ Management gas flow: Max. 65 kPa ▶ Applicable to hospital pipeline CO2 ▶ Change the flow level on the front panel (Standard and low flow) ▶ 6 kg ▶ 145 mm (W) x 170 mm (H) x 390 mm (D)

Cleaning Brushes

Endoscopes		Where to use	Suitable brush
ELUXEO™ 700 Series	Reusable	Working channel	WB7024FW
	Reusable	Cylinder port	WB11003FW
600 Series	Reusable	Working channel	WB4321FW2
	Reusable	Cylinder port	WB11002FW2
500 Series	Reusable	Working channel	WB4321FW2; WB3530FW
	Reusable	Cylinder port	WB11002FW2
Duodenoscopes			
ED-530XT/XT8	Reusable	Working channel	WB4321FW2
	Reusable	Cylinder port	WB11002FW2
	Reusable	Air-water channel	WB2221FW2
	Disposable	Working channel	WB7025DC
	Disposable	Cylinder port	WB11003DV
	Disposable	Elevator mechanism	WB1318DE
Short Double-Ballon Endoscope			
EI-580BT	Reusable	Working channel	WB7024FW
	Reusable	Cylinder port	WB11003FW

Endoscopes		Where to use	Suitable brush
Bronchoscopes			
EB-580	Reusable	Working channel	WB7024FW
	Reusable	Cylinder port	WB11003FW
EB-530 series	Reusable	Working channel	WB3212FW2; WB3503FW
	Reusable	Cylinder port	WB11002FW2
Ultrasonic Videoscopes			
EG-580UT/UR	Reusable	Working channel	WB7024FW
	Reusable	Cylinder port	WB11003FW
	Reusable	Balloon channel	WB2221FW2
	Disposable	Elevator mechanism	WB1318DE
EG-530UT2/UR2	Reusable	Working channel	WB4325DS
	Reusable	Cylinder port	WB11002FW2
EG-530UR2	Reusable	Working channel	WB3225DS
EB-530US	Reusable	Working channel	WB3212FW2; WB3503FW
	Reusable	Cylinder port	WB11002FW2
Fibrescopes			
FG-1Z	Reusable	Working channel	WB3518FW
FC-1Z	Reusable	Working channel	WB3524FW

MONITORS

Model	Product Description
-------	---------------------

Monitors*	
Radiance 27" ULTRA 	 Image: Radiance 27" Including FUJIFILM BIOS for the best performance
Radiance 26" G2 HB LED 	
Radiance 24" G2 LED 	
Radiance 19" HD 	
	27" HD type LCD monitor with Ultra bright LED Backlight ▶ High-Definition, Color Correction Technology (CCT), Full Multi-Modality, Gorilla Glass front panel ▶ 8.9 kg, 678 mm (W) x 445 mm (H) x 84 mm (D) ▶ Input signal: HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA ▶ Output signal: HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
	26" HD type High Brightness LCD monitor with LED Backlight ▶ High-Definition, Color Correction Technology (CCT), Full Multi-Modality ▶ 8.2 kg, 673 mm (W) x 418 mm (H) x 88 mm (D) ▶ Input signal: HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA ▶ Output signal: HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
	24" HD type LCD monitor with LED Backlight ▶ High-Definition, Color Correction Technology (CCT), Full Multi-Modality ▶ 7.1 kg, 597 mm (W) x 401 mm (H) x 100 mm (D) ▶ Input signal: HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, Composite, VGA ▶ Output signal: HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
	19" HD type LCD monitor ▶ High-Definition, Full Multi-Modality ▶ 6.8 kg, 465 mm (W) x 400 mm (H) x 98 mm (D) ▶ Input signal: HD-SDI x 2, DVI-D, DVI-I, RGBS, YPbPr, S-Video, composite, VGA ▶ Output signal: HD-SDI, DVI, RGBS, YPbPr/VGA, S-Video, Composite
EndoVue 19" 	 19" HD type LCD monitor with LED Backlight ▶ 4,2 kg, 464.8 mm (W) x 396.2 mm (H) x 99 mm (D) ▶ Input Signal: DVI-D, HD-SDI, (HD-)RGBS, (HD-)YPbPr, VGA, S-Video, Composite


*Manufactured by NDS Surgical Imaging, LLC, USA



CARTS

Model	Product Description
Carts*	
ITD-03	 Compact cart with trafa ▶ Approx. 76 kg, 618 mm (W) x 1,331 mm (H) x 624 mm (D) ▶ 1x endoscope hanger, 3x shelves, 1x monitor stand ▶ 1,200 VA, Isolating transformer Image: Cart ITD-03 plus
ITD-03 GAS	 Uni-cart CO2 insufflator ▶ 35 kg, 540 mm (W) x 580 mm (H) x 857 mm (D) ▶ 3x shelf, 1x gas cylinder holder (Ø145 mm, max. 20 kg) Image: Cart ITD-03 GAS
ITD-03 Sono	 Uni-cart sonoprobe ▶ 40 kg, 540 mm (W) x 1,510 mm (H) x 580 mm (D) ▶ 3x shelves, 1x monitor stand ▶ 1x jointed arm ▶ 660 VA, Isolating transformer Image: Cart ITD-03 Sono

*Manufactured by ITD GmbH, Germany

Model	Product Description
Carts*	
ITD-03 Iso	 Compact cart with insulation monitoring function ▶ Approx. 76 kg, 618 mm (W) x 1,331 mm (H) x 624 mm (D) ▶ 1x endoscope hanger, 3x shelves, 1x monitor stand ▶ 2,000 VA, Isolating transformer with Insulation monitoring function Image: Cart ITD-03 Iso plus (with insulation monitoring function)
ITD-03 EUS	 Compact cart EUS system ▶ Approx. 85 kg, 618 mm (W) x 1,491 mm (H) x 624 mm (D) ▶ 2x endoscope hanger, 4x shelves, 2x monitor stand ▶ 2,000 VA, Isolating transformer with Insulation monitoring function Image: Cart ITD-03 EUS Iso
ITD-03 EUS with upgrade Kit	 Compact cart EUS system with upgrade Kit ▶ Approx. 88 kg, 818 mm (W) x 1,765 mm (H) x 624 mm (D) ▶ 2x endoscope hanger, 4x shelves, 2x monitor stand ▶ 1x side shelf with jointed arm ▶ 2,000 VA, Isolating transformer with Insulation monitoring function Image: Cart ITD-03 EUS Iso
ITD-00+ Also available without back door	 Classic cart ▶ Approx. 70 kg, 720 mm (W) x 1630 mm (H) x 720 mm (D) ▶ 1x endoscope hanger with one-hand-adjustment, 1 shelf with a handle, 1 monitor stand Image: Cart ITD-00 with back door



FUJIFILM

FUJIFILM Europe GmbH

Heesenstr. 31, 40549 Düsseldorf, Germany
Tel.: +49 211-50 89 0, Fax: +49 211-50 89 8700
www.fujifilm.eu, endoscopy_eu@fujifilm.com

Specifications are subject to change without notice. The name FUJIFILM and the FUJIFILM logo are trademarks of FUJIFILM Corporation. All other trademarks shown are trademarks of their respective owners. All rights reserved. 07/2018