

# GW-100

Endoscopic CO<sub>2</sub> Insufflator

NEW



FUJIFILM


- ▶ Fast resorption of insufflated CO<sub>2</sub> for timesaving and patient friendly examinations

Our latest GW-100 CO<sub>2</sub> insufflator offers clinicians an optimized and easy to handle procedure as well as a maximum of patients comfort.

## Features

- ▶ Direct connection to hospitals medical CO<sub>2</sub> pipeline as well as to medical CO<sub>2</sub> cylinder
- ▶ Easy-to-use CO<sub>2</sub> flow rate switching function and compact design
- ▶ 2 controlled flowrate settings

## System configuration



Tube sets for the connection of GW-100 to the medical gas pipeline and medical gas cylinders are available.

## Specifications

- ▶ GW-100

Dimensions	145 x 170 x 390 mm (W x H x D)
Weight	6 kg
Power Requirements	100 - 240VAC ± 10%, 50/60Hz, 0.3A
Applicable Gas	CO <sub>2</sub> gas for medical use (Connectable to CO <sub>2</sub> cylinder or medical gas pipeline)
Gas Supply Pressure	65kPa (max.)

